

deSoutter
MEDICAL



Accessories

FOR SURGICAL POWERED INSTRUMENTS

www.de-soutter.com

CONTENTS

Supplied Accessories

Chuck keys	3
Wrenches	3
Lubricating oil	3
Cleaning brush	3

Battery and Electrical Components

Aseptic batteries	4
Sterile batteries	4
Battery charger.....	5
Power supply units	5
Foot pedals	5

Sterilisation equipment

Sterilisation cases and insert trays.....	8
---	---

Pneumatic instrument accessories

“Wash Assist” air hoses	11
Standard air hoses	12
Air hose parts	13
Extension hoses	14
Conversion hoses	14
Air fittings	14
Multi stage regulator	14
Air foot control	15
Dermatome	15
Mains supply cable	15








General accessories

Chucks and adaptors	16
---------------------------	----

Sales and Service Policies

Servicing and repairs.....	18
Warranty	19

SUPPLIED ACCESSORIES

CHUCK KEYS, WRENCHES AND SPANNERS			
	Used with	Description	Order code
	Q series adaptors, DX-600/500, DU-470, DZ-450/400, DV-225/220/200	Key for 0.5 - 4 mm chuck	8780
	Q series adaptors, RX-600/500, RU-470, DX-600/500, DZ-450/400, RZ-450/400, DU-470	Key for 0.5 - 7.4 mm chuck	30062
	D32-100, MV-220	Sagittal saw blade screwdriver	2250
	D34-100	Reciprocating wrench	2260
	CU-470, CZ-450/400, CV-220/200, BV-220/200	Wrench	9940
	ECV-221/200, BU-470, BZ-450, BV-220	Wrench	14770
	SV-220/200, MV-220/200, GD-113	Saw blade hexagon driver	9930
	SZ-450/400, SU-470	Saw blade hexagon driver	9590
	DPX-170, NS3	12 mm/22 mm spanner	190723
WIRE GUARDS			
	D36-100	Wire guard	10360
	For use with small bone (MCZ, MBU) wire and pin drivers	Wire guard	9570
LUBRICANTS			
	For use with Large and Small bone instruments - MPX, MPZ, DPX, EPV.	25 ml bottle of mineral oil	30982
	For use with Large and Small bone instruments - DPZ, D30, PLX.	50 ml bottle of mineral oil	12570
BRUSHES			
	WZ-451, WU-470/471, PU-470, PQ-701/702, WQ-701/707, PZ-450	Cleaning brush	8120
	WX-630, WQ-700/710	Cleaning brush	15090
















BATTERY AND ELECTRICAL COMPONENTS

For use with powered handpiece types MBQ, DBK, MBU, DBN, DBC, DBR, MDX and TDZ

STERILE BATTERIES

	SB-703 Small Li-Ion sterile battery 13.2V Capacity 1100 mAh Order code 17210		SB-704 Large Li-Ion sterile battery 13.2V Capacity 2500 mAh Order code 17220		SB-600 9.6 Vdc NiCd Capacity 1200 mAh Battery Pack Order code 11630
---	---	---	---	---	---

ASEPTIC BATTERIES

	AB-703 Li-Ion Small aseptic battery 10.95V Capacity 2000 mAh Order code 17910		AH-703 Li-Ion Small aseptic battery housing Order code 18080		AS-703 Li-Ion Small aseptic battery shield Order code 18100
	AB-704 Li-Ion Large aseptic battery 14.6V Capacity 2000 mAh Order code 17920		AH-704 Li-Ion Large aseptic battery housing Order code 18090		AS-704 Li-Ion Large aseptic battery shield Order code 18110
	AB-450 Small aseptic battery 10.8V NiMH Capacity 1000 mAh Order code 14510		AH-450 Small aseptic battery housing Order code 14520		AS-450 Small aseptic battery shield Order code 14530
	AB-450 Small aseptic battery 10.8V NiMH Capacity 1000 mAh Order code 14510		AH-700 Small aseptic battery housing Order code 16140		AS-700 Small aseptic battery shield Order code 16150
	AB-600 Large aseptic battery 9.6V NiMH Capacity 2100 mAh Order code 11520		AH-610 Large aseptic battery housing Order code 14430		AS-610 Large aseptic battery shield Order code 14420
	AB-701 Medium aseptic battery 10.8V NiMH Capacity 1950 mAh Order code 15340		AH-701 Medium aseptic battery housing Order code 15670		AS-701 Medium aseptic battery shield Order code 15890
	AB-702 Large aseptic battery 14.4V NiMH Capacity 1950 mAh Order code 15350		AH-702 Large aseptic battery housing Order code 15680		AS-702 Large aseptic battery shield Order code 15900

Battery charger, power supply, power consoles and foot pedals





BATTERY CHARGERS



BC-700 Multi bay charger

* CA-450 charging adaptor required
Order code 14540



Model	Order code	Battery Type	Charger Module	Module Order code
BC-700	15700	-	None	-
	16200	AB-450*/AB-600 SB-600	CA-702	15770 
	16210	AB-701/AB-702	CA-701	15760 
	16220	SB-703/SB-704	CA-700	15750 
	16230	AB-703/AB-704	CA-703	18120 
BC-705	17980	AB-701/AB-702	-	-
BC-706	17990	SB-703/SB-704	-	-
BC-707	18000	AB-703	-	-
BC-708	18010	AB-704	-	-



Single bay charger

POWER SUPPLY



PS-701
Corded direct power supply.
Used with MBQ and MBU.
Power supply input 240 Vac.
Order code 17870



CM-700
Corded module
Order code 15720

POWER CONSOLES



PC-470
Used with MCI, MCU
handpieces
Order code 18050



PC-471
less irrigation
Used with MCI, MCU
handpieces
Order code 19080



PC-120 less irrigation
Order code 103305
PC-121 with irrigation
Order code 103318
Used with HDS-121



PC-472
with irrigation
Used with MCI, MCU
handpieces
Order code 19090


4 metre cable for
connecting
handpiece to
PC-470, PC-471,
PC-472
Power Console
Order code 18570

Power supply input 230Vac, 120Vac and 100Vac.


FOOT PEDALS



FS-270
Order code 18580
Used with
MCI, MCU



Wired
MFP-120
Order code 103556
Used with
HDS-121




Wireless
MFP-121
Order code 103562
Used with
HDS-121

STERILISATION EQUIPMENT


Sterilisation containers

STANDARD SIZE WIRE CONTAINER



Standard wire container - Order code 19010
Accepts standard insert tray

Standard wire insert tray
Order code 668033



Silicone mat for tray -
Order code 616153

HALF SIZE WIRE CONTAINER



Half size wire container
Order code 19000
Accepts half size insert tray

Half size wire insert tray
Order code 667923



STERILE FILTERED CONTAINER



Standard Sterile filtered container
Order code 19250

PTFE sterile filtered container
Order code 19810

Accepts standard wire
insert tray

Lid gasket
Order code 670023



PTFE filter cassette
Order code 671073

Filter paper - single use
pack of 100
Order code 657243



Security tag,
pack of 100
Order code 657253

Sterile indicator label, pack of 100 Order code 657263

HALF SIZE STERILE FILTERED CONTAINER




Half size sterile Filtered Container
Order code 19370

Accepts half size wire
insert tray

Lid gasket -
Order code 669903


Large bone battery sterilisation systems

ORTHODRIVE® BATTERY SYSTEM




Set up	System	Wire Container and insert tray Order code	Insert tray Order code
3 MBQ/DBK/DBC handpieces 9 attachments 2 aseptic housings and shields	Aseptic Li-Ion	18150	17960
2 MBQ/DBK/DBC handpieces 9 attachments 2 aseptic housings and shields	Aseptic Li-Ion	18160	17970
3 MBQ/DBK/DBC handpieces 9 attachments, 2 batteries	Sterile	16250	17110
2 MBQ/DBK/DBC handpieces 9 attachments, 2 batteries	Sterile	19140	19040

ORTHODRIVE® LITE BATTERY SYSTEM



Set up	System	Wire Container and insert tray Order code	Insert tray Order code
1 MBQ-707/708 handpiece 8 attachments, 8 adaptors 1 aseptic housings and shield	Aseptic Li-Ion	18130	17900
1 MBQ-707/708 handpiece 8 attachments, 8 adaptors 1 sterile battery	Sterile	17830	17550
2 MBQ-707/708 handpieces 6 attachments, 8 adaptors 2 aseptic housings and shields	Aseptic Li-Ion	19120	19020
2 MBQ-707/708 handpieces 6 attachments, 8 adaptors 2 sterile batteries	Sterile	19160	19060
1 MBQ-707 handpiece 4 attachments, 1 aseptic battery, housing and shield	Aseptic Li-Ion	19400 half size	19290 half size

MULTIDRIVE® PNEUMATIC SYSTEM



Set up	System	Wire Container and insert tray Order code	Insert tray Order code
2 MPX/DPX handpieces	Pneumatic	14640	14650




Washing cap for MPX, SPX and DPX - Order code 9950

Insert trays shown complete with handpieces and other items are for illustration purposes only.

STERILISATION EQUIPMENT


Trauma and Small Bone sterilisation systems

TRAUMADRIVE™ BATTERY SYSTEM




Set up	System	Wire container and insert tray Order code	Insert tray Order code
2 MBU handpieces 9 attachments 2 aseptic housings and shields	Aseptic Li ion	18140	18060
1 MBU handpiece 9 attachments 2 aseptic housings and shields	Aseptic Li ion	19130	19030
2 MBU handpieces 9 attachments 2 sterile batteries	Sterile	16940	16980
1 MBU handpiece 9 attachments 1 sterile battery	Aseptic Li ion	19150	19050
1 MBU handpiece 6 attachments 1 aseptic housing and shield	Aseptic Li ion	19210 Half size	19190 Half size

OSTEODRIVE™ ELECTRIC SYSTEM



Set up	Wire Container and insert tray Order code	Insert tray Order code
1 MCI-270 handpiece 1 MCU-470 handpiece 6 MCI-270 attachments 9 MCU-470 attachments 2 connection cables	18950	18990

MULTIDRIVE™ PNEUMATIC SYSTEM




Set up	System	Wire Container and insert tray Order code	Insert tray Order code
1 MPZ and 1 EPV handpiece, 6 EPV attachments, 9 MPZ attachments and air hose	Pneumatic	14390	14410



Washing cap for EPV, MPZ handpieces - Order code 9950

MICRODRIVE™ - PNEUMATIC SYSTEM



D30 Series Handpieces
Order code 7540




Washing cap for D30 Series and PLX handpieces
Order code 8110

Insert trays shown complete with handpieces and other items are for illustration purposes only.


Cardiothoracic and battery systems

STERNUDRIVE™ ASEPTIC BATTERY SYSTEM




Set up	Half size wire container and insert tray Order code	Insert tray only Order code
1 DBR-700 handpiece 1 aseptic housing 1 aseptic shield	19170	19070
1 DBR- 470/DBC-470 handpiece 1 sternum guard (DBC-470 only) 1 aseptic housing and shield	18210	18200
1 MBC-709 handpiece 3 attachments 1 sternum guard 1 aseptic housing and shield	19280	19300

STERNUDRIVE™ STERILE BATTERY SYSTEM



Set up	Half size wire container and insert tray Order code	Insert tray only Order code
1 DBR-700 handpiece 1 sterile battery	17890	17880
1 DBR- 470/DBC-470 handpiece 1 sternum guard (DBC-470 only) 1 sterile battery	19270	19110
1 MBC-709 handpiece 3 attachments 1 sternum guard 1 sterile battery	19200	19100

STERNUDRIVE™ PNEUMATIC SYSTEM



Set up	Half size wire container and insert tray - Order code	Insert tray only Order code
1 DPC-470 handpiece 1 sternum guard, 1 hose	17160	17170
DPX-170	Standard size wire container, lid and tray - Order code 19010 Silicone mat for tray - Order code 616153	




Washing cap for MPX and DPX Order code 608073



Washing cap for DPC-470 Order code 9950

BATTERY INSERT TRAY




Set up	Half size wire container and insert tray Order code	Insert tray only Order code
2 off SB-703 or SB-704 sterile batteries	17390	17400

Insert trays shown complete with handpieces and other items are for illustration purposes only.

STERILISATION EQUIPMENT

Neuro, ENT and Dermatome systems

HEKADRILL™ INSTRUMENT SYSTEM



Set up - electric instrument tray and bur head tray	Half size wire container, electric tray and bur head tray, Order code
1 HDS-121 handpiece 1 P-P01 perforator 1 electric cable 5 attachments and spanner	103610
Set up - pneumatic instrument tray and bur head tray	Half size wire container, pneumatic tray and bur head tray, Order code
1 HDS-120 handpiece 1 P-P02 perforator 1 pneumatic hose 5 attachments and spanner	103629


Tray only for electric system - Order code 103953
 Tray only for pneumatic system - Order code 103989
 Tray only for bur heads - Order code 103966
 Wire sterilisation container and lid - Order code 103944

 Washing cap for HDS-120 and P-202 handpieces. Order code 9950





Sterile filtered container
 Accepts half size electric or pneumatic instrument trays plus bur head insert tray
 Order code 103977


NEURO/ENT PNEUMATIC SYSTEM



Model	Container and insert tray Order code
HST/DPZ	12550

Plastic base, plastic lid and plastic insert
 Washing cap for HST/DPZ handpieces. Order code 9950
 Washing cap for HST nose only. Order code 9950

DERMADRIVE™ ELECTRIC SYSTEM



Model	Container and insert tray Order code
GD-113 plus accessories	8570

Plastic base, plastic lid and plastic insert











PNEUMATIC INSTRUMENT ACCESSORIES

Wash Assist air hoses

WASH ASSIST AIR HOSES



The Wash Assist co-axial hose system is recommended where machine washing of air hoses is demanded. By releasing the probe end from hose, which is retained on the hose via the rubber strap, the hose ends can be joined and sealed prior to washing.

Probe end	For model type	Complete hose length - Order code		
		3 metre	4 metre	5 metre
 MA7-Schrader	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13000	13100	13200
		13300	13400	13500
		13600	13700	13800
 A.O.	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13010	13110	13210
		13310	13410	13510
		13610	13710	13810
 AGA Inst.	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13020	13120	13220
		13320	13420	13520
		13620	13720	13820
 Draeger	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13030	13130	13230
		13330	13430	13530
		13630	13730	13830
 Heyer	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13040	13140	13240
		13340	13440	13540
		13640	13740	13840
 Air Liquide	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13050	13150	13250
		13350	13450	13550
		13650	13750	13850
 DIN	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13060	13160	13260
		13360	13460	13560
		13660	13760	13860
 AGA Air	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13080	13180	13280
		13380	13480	13580
		13680	13780	13880
 Flowtech	MPX, DPX, SPX, MPZ, EPV, HST,, DPC-470 PLX, D30 Series	13070	13170	13270
		13370	13470	13570
		13670	13770	13870
 NP05	MPX, DPX, SPX, MPZ, EPV, HST, DPC-470 PLX, D30 Series	13090	13190	13290
		13390	13490	13590
		13690	13790	13890

To special order






PNEUMATIC INSTRUMENT ACCESSORIES


Standard air hoses and wash caps

STANDARD AIR HOSES



Quick release, bayonet action hose fitting with combined coaxial air supply and exhaust hose. This system is non-preferred for machine washing.

Probe end	For model type	Complete hose length - Order code		
		3 metre	4 metre	5 metre
 MA7-Schrader	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470 PLX, D30 Series	5300 10900 7500	5310 10910 9320	5320 10920 7640
 A.O.	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470	5330 10930	5340 10940	5350 10950
 AGA Inst.	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470 PLX, D30 Series	5360 10960 8720	5370 10970 8730	5380 10980
 Draeger	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470 PLX, D30 Series	5390 10990 8580	5400 11000 8590	5410 11010 8600
 Heyer	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470 PLX, D30 Series	7820 11020 8750	7830 11030	7840 11040 8770
Air Liquide	MPX, DPX, SPX, HDS-120 MPZ, EPV, HST, DPC-470	10390 11050	10400 11060	10410 11070

 To special order

HOSE WASHING CAP



Probe end	For model type		
	MPX, DPX, SPX, Order code	MPZ, EPV, HST, DPC-470 Order code	PLX, D30 Series Order code
MA7 Schrader	8950	9650	11430
A.O.	10880	9670	11450
AGA		9660	
Draeger	10890	9680	11460
Heyer		9660	

Wash Assist air hose parts, adaptor without air hose, air hose without adaptor

WASH ASSIST AIR HOSE PARTS



Wash assist co-axial hose.

Hose ends connected with probe detached and retained with rubber strap.



Replacement rubber strap

Order code 632393

AIR HOSE WITHOUT ADAPTOR



For model type	Hose length - Order code		
	3 metre	4 metre	5 metre
MPX, DPX, SPX,	632303	632313	632323
MPZ, EPV, HDS-120, P-P02, HST	632333	632343	632353
DPC-470			
PLX, D30 Series	632363	632373	632383

ADAPTOR WITHOUT HOSE




Adaptor type	For model type - Order code		
	MPX	MPZ	D30
MA7	632003	632103	632203
A.O.	632013	632113	632213
AGA Inst.	632023	632123	632223
Draeger	632033	632133	632233
Heyer	632043	632143	632243
Air Liquide	632053	632153	632253
DIN	632063	632163	632263
Flowtech	632073	632173	632273
AGA Air	632083	632183	632283
NPO5	632093	632193	632293

PNEUMATIC INSTRUMENT ACCESSORIES

Extension hoses, conversion hoses, convertors and air regulator.


EXTENSION HOSES



Hose length - Order code	
3 metre	6 metre
5670	6120


Hose length - Order code	
3 metre	5 metre
8940	10370

CONVERSION HOSES



Probe end	Female end	Order code
A.O.	MA7	5150
AGA Inst.	MA7	7510
Draeger	MA7	7520
Air Liquide	MA7	8840
1/4" NPT	MA7	6290

CONVERTORS AND AIR REGULATOR



MA7 Probe to double MA7 female
Order code 10420



MA7 Probe to A.O. female
Order code 12400

A.O. Probe to MA7 female
Order code 12410

MA7 female to 1/4" NPT female thread
Order code 5160

1/4" BSPT male to 5/8" x 27 TPI male thread
Order code 5170



Multi stage regulator with MA7 female outlet - bull nose

Inlet pressure	Outlet pressure	Air Inlet	Order code
21 - 200 bar	0 - 10 bar	5/8" BSP	5180

PNEUMATIC INSTRUMENT ACCESSORIES AND DERMATOME

Air foot control, Dermatome and mains cable

AIR FOOT PEDAL - for use with EPV-221,MPZ-451,HST-100,HDS-120,DPZ,P-P02



Model	Order code	Foot pedal features	
		Inlet connection	Outlet connection
FV-220	14900	Integral 3 metre supply hose with MA7 probe	Accepts MPZ air hoses - see page 14
FV-221	14910	Accepts Draeger Wash Assist hoses with air being exhausted at wall fitting	Accepts MPZ air hoses - see page 14
FV-222	15000	Integral 3 metre supply hose with MA7 probe	MA7 outlet

DERMATOME



Model	230 V UK plug Order code	230 V Euro plug Order code	120 V US plug Order code
CM2 motor controller, micro motor and foot pedal	1292334	1292334	1292334



Sole plate 23 mm cutting width Order code 17340
 Sole plate 48 mm cutting width Order code 17260
 Sole plate 73 mm cutting width Order code 17270

MAINS CABLE



Used with:
 Power Console
 BC-700 Battery Charger
 PS-701 Power Supply
 IPV-220 Irrigation Pump

















Cable length 3 metres

UK	Order code 8290
Euro	Order code 8300
US/Canada	Order code 8310
Australia	Order code 8320
New Zealand	Order code 8320
Swiss	Order code 14700

GENERAL ACCESSORIES

All drives on this page are cannulated excluding A.O. drill types marked *

CHUCKS AND ADAPTORS

Input drive										
Output drive		Q Series	Hudson	Zimmer	Harris	Jacobs	DIN 58809	A.O. Reamer	Trinkle	A.O. Drill
	Jacobs 4mm	15360								6270*
	Jacobs 7.4mm (1/4")	15370	5440	8800	8810		8820	7970	8830	11550*
	3 mm (1/8") Keyless	15380	12320	12430				12310		
	7.4 mm (1/4") Keyless	15390	12340	12440				12330		
	ISO 6360									11270*
	A.O. Drill	15480	5470	10260	10270	10280	10290	7980	10300	
	A.O. Alesatore	15490	5480	10170	10180	10190	10200		10210	
	Hudson	11530		10000	10010	10020	10030	6560	10040	11580*
	Zimmer	15460	5450		10050	10060	10070	6590	10080	
	Trinkle	15540	7910	10220	10230	10240	10250	6580		
	Harris	17130	5460	10090		10100	10110	6550	10120	
	DHS	15500	11260					7770		
	DCS		17510					17930		
	Targon TX	17930								
	DIN 58809		5490	10130	10140	10150		6570	10160	
	Biomet drill									12890*



SALES AND SERVICE POLICIES

Servicing and Repairs

Our equipment is supported by a number of strategically located servicing and repair centres. All powered instruments and attachments should be regularly serviced to ensure optimum performance and reliability. Regular servicing will minimise the risk of unexpected breakdowns. Servicing includes a comprehensive functional check, cleaning and testing, followed by the replacement of any worn or damaged components. Because of specialist techniques used in the manufacture and maintenance of our equipment, no user or third party servicing is possible.

To return equipment for repair or warranty, please ensure the following:

- Disinfect/sterilise the equipment in accordance with reprocessing instructions outlined in the products Operating Instructions.
- Record the serial number of the equipment being returned with the return documentation and enclose a brief statement describing the reason for the equipment's return together with a contact information. (contact name, email address and telephone number).
- For warranty returns, please state the original invoice number, date of purchase and a brief statement describing the reasons for the equipment's return.
- For repair/servicing, please enclose a purchase order number to cover the Minimum Repair/Inspection Charge.
- Pack instruments securely, as any in-transit damage is not covered by the warranty.
- The cost of all repairs will always be advised by De Soutter Medical and approved by the customer before any work commences.

Lithium Ion batteries - returns and transportation

Please see individual Operating Instructions for our policy for disposal, return and transportation of Lithium Ion batteries. Lithium Ion battery packs should not be returned or transported if there are any signs of damage. If in doubt, please seek guidance from your local De Soutter Medical representative.

- When returning a battery pack, our approved Lithium Ion battery label should be printed from the De Soutter Medical website and attached to the outer packaging. For more information on Lithium Ion battery transportation, please visit our support page "www.de-soutter.com/Support"

Returned goods

Requests for credit returns must be received within 60 days of the invoice date.

Equipment returned to De Soutter Medical for credit without our prior written authorisation will not be accepted. All returned goods will be subject to a re-stocking charge of 15% of the net original invoice value.

Credit cannot be given for items damaged during return shipment due to inadequate packaging. To receive a full credit, less the re-stocking charge, all products must be returned in full package quantities and any consumables must be in the original, undamaged packaging.

Warranty

De Soutter Medical, (“the Company”), warrants to the first purchaser or lessee (“Customer”) that the instruments, attachments, parts and consumables manufactured by or for the Company (hereinafter collectively called “Products”) are guaranteed to be free from defects in material and workmanship for the warranty period for a particular Product, beginning from the date of invoice.

New Products - Warranty Period

- Powered Surgical Handpieces and Attachments - Twelve (12) months.
- Cast and Autopsy Saws - Twelve (12) months.
- Battery Chargers and Power Supply/Control Units - Twelve (12) months.
- Pneumatic Hoses - Twelve (12) months.
- Handpiece cords and Power Cords - Twelve (12) months.
- Electric and Pneumatic Footswitches - Twelve (12) months.
- Handpiece Accessories, Bur Guards, Blade Guards and Adaptors - Twelve (12) months.
- Battery Packs - Six (6) months.
- Consumables i.e. single use Burs, Blades and Pump Kits - for such period expiring on the “use by date”.

This limited warranty is restricted to repair or replacement by the Company, at its option, of any Product found to be defective during the warranty period. All warranties shall apply only to the original Customer. In no event shall the Company be liable for any anticipated profits, consequential damages or loss of time incurred by the Customer with the purchase or use of any Product.

The Company shall be under no liability in respect of any defect arising or resulting from fair wear and tear, improper or inadequate installation, use or maintenance, accidental or wilful damage, negligence, abnormal working conditions, failure to follow the Company’s product operating instructions, misuse or modification or alteration or repair of the Products by unauthorised third parties or the Customer.

Repaired Products

The Company warrants to the Customer that any service or repair work performed by the Company on Products not within the warranty period shall be free from defects in workmanship and materials for a period of twelve (12) months, six (6) months for pneumatic hoses, after the date the Customer received the Product(s).

The aforementioned warranty does not affect your Statutory Rights.



UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited
Halton Brook Business Park, Weston Road
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Fax: +44 (0) 1296 634 033
Email: info@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Fax: +61 (0) 3 9702 4484
Email: australia@de-soutter.com

AUSTRIA

De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Fax: +43 (0) 2243 21 656
Email: austria@de-soutter.com

BELGIUM

De Soutter Medical Belgium
Bessemmerstraat 14
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Fax: +32 (0) 89/70 12 19
Email: belgium@de-soutter.com

FRANCE

De Soutter Medical France
1 252 Avenue du Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Fax: +33 (0) 9 70 61 37 60
Email: france@de-soutter.com

GERMANY

De Soutter Medical Deutschland
Bahnhofstraße 4
D-66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Fax: +49 (0) 68 52-99 12 47
Email: deutschland@de-soutter.com

ITALY

De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Fax: +39 (0) 2 9009 2673
Email: italy@de-soutter.com

NETHERLANDS

De Soutter Medical Nederland
Gelderlandhaven 2X
3433 PG Nieuwegein
Tel: +31 (0) 85 0491480
Fax: +31 (0) 85 0491489
Email: nederland@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
Mooresville, NC 28117
Tel: +1 (704) 655 9040
fax: +1 (704) 987 2035
Email: usa@de-soutter.com

deSoutter
MEDICAL

For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.